Disinfection cold fog device ST-83.1

The brand new ST-83.1 cold fog hand tool is ideal for small applications for fogging of disinfectants, odor neutralizers or insecticides.

- » Inexpensive and safe application for hygienically germ-free cleaning of interiors to remove microorganisms, viruses and bacteria. Also for use on construction sites for dust deposition or mold removal
- » No inhalation of the disinfectant solution as rooms do not have to be entered
- » The finest nebulization of active ingredients comes into the smallest gaps and into the last corner
- » A compressed air connection is required: Inlet pressure 3 bar, 180-220 I/min.
- » Extremely productive use: (e.g. 1 liter for rooms with 60 m³, room size of approx. 25 m² or 5 - 7 vehicle compartments



Disinfection cold fog air device ST-83. For compressed air connection 3 bar, 180-220 l / min. Stainless steel / plastic

R+M Nr.	D		
200 083 500	0:7	1 l	









Two different fog nozzle inserts are included.



Particularly suitable in smaller rooms such as passenger compartments.





Larger throw range. For the treatment of entire rooms or for targeted treatment at the air intake of ventilation systems.



ST-83 cold fog hand tool accessories

ST-83 window bracket



With elastomer protection rings. Stainless steel

Starriess steer	
R+M Nr.	
200 083 400	

ST-83 adapter for canisters



For 5 - 30 liter canisters. Incl. 0.5 m clear PVC hose Plastic. 10 nozzle inserts (0.5 - 2.0 mm). for suction of chemicals, DN 4

Tot suction of enermedis, bit i				
	R+M Nr.			
	200 083 200			

Exchangeable nozzle inserts



Incl. o-rinas

3
R+M Nr.
200 163 340